NAPA-SOLANO-YOLO-MARIN COUNTY PUBLIC HEALTH LABORATORY

2201 COURAGE DRIVE, MS 9-200 FAIRFIELD, CA 94533 (707) 784-4410 FAX (707) 423-1979



Test	Sarcoptes scabiei (scabies) screening
Test Description	Detection of Sarcoptes scabiei uses microscopic analysis of skin scarping to detect the potential parasitic mite via saline mount.
Acceptable specimens	Skin scrapings (preferably between fingers and fold of wrist) should be collected in a dry, sterile test-tube or in a clean paper wrap.
Specimen storage and shipping	Keep the specimen at room temperature and transport to the laboratory as soon as possible.
CPT code	
Test fee	
Result availability (turnaround time)	24hr